

**Worksheet for the determination of the absorbed dose to water  
in an electron-beam**

User: **CLINICA AUNA CHICLAYO** Date: **31-Mar-25**

**1. Radiation treatment unit and reference conditions for  $D_{w,Q}$  determination**

Accelerator:	<b>INFINITY CHICLAYO</b>			
Nominal dose rate:	<b>600.0</b>	MU min <sup>-1</sup>	Nominal energy:	<b>12.000</b> MeV
Reference phantom:			Measured $R_{50}$ :	<b>4.740</b> g cm <sup>-2</sup>
Reference field size:	<b>10X10</b>	cm x cm	obtained from	
Beam quality, $Q (R_{50,w})$	<b>4.817</b>	g cm <sup>-2</sup>	Reference SSD:	<b>100</b> cm
			Ref. depth $z_{ref,w} = 0.6 R_{50} - 0.1$ :	<b>2.8</b> g cm <sup>-2</sup>

**2. Ionization chamber and electrometer**

Ion. chamber model		Serial No.:	<b>268895004</b>
Chamber wall / window material:	<b>PMMA</b>	thickness:	<b>0.1180</b> g cm <sup>-2</sup>
Waterproof sleeve material:		thickness:	
Phantom window material:		thickness:	

**Abs. dose-to-water calibration factor <sup>a</sup>**

Calibration quality $Q_0$ :		Calibration depth:	<b>5.0</b> g cm <sup>-2</sup>
If $Q_0$ is electron beam, give $R_{50,w}$ :			
Reference conditions for calibration			
$P_0$ :	<b>101.3</b> kPa	$T_0$ :	<b>20.0</b> °C
		Rel. humidity:	<b>50</b> %

Pol. potential $V_1$ :		V	Calib. polarity:	
			User polarity:	

Calibration laboratory:	<b>LSCD IPEN</b>	Date:	<b>26-Jun-24</b>
Electrometer model:	<b>PC ELECTROMETER</b>	Serial no.:	<b>270267006</b>
Calib. separately from chamber:		Range setting:	
If yes Calibration laboratory:		Date:	

**3. Phantom**

Water phantom window material:		thickness:		g cm <sup>-2</sup>
Plastic phantom phantom material:		density:		g cm <sup>-3</sup>
depth scaling factor $c_{pl}$ :		reference depth $z_{ref,pl} = z_{ref} / c_{pl}$ :		g cm <sup>-2</sup>
fluence scaling factor :		$h_{pl} =$		

**4. Dosimetry reading <sup>b</sup> and correction for influence quantities**

Uncorrected dosimeter reading at $V_1$ and user polarity:	<b>12.788</b>
Corresponding accelerator monitor units:	<b>100</b> MU
Ratio of dosimeter reading and monitor units:	$M_1 =$ <b>0.1279</b>
(i) $P$ :	<b>100.6</b> kPa
$T$ :	<b>23.5</b> °C
Rel. humidity:	<b>50</b> %
	<b>1.0189</b>

(ii) Electrometer calibration factor  $k_{elec} =$

(iii) Polarity correction <sup>c</sup>	rdg at $+V_1$	$M_+ =$ <b>12.788</b>	rdg at $-V_1$ :	$M_- =$ <b>-12.794</b>
--	---------------	-----------------------	-----------------	------------------------

**1.0002**

  
 Clinicas **auna** | Sede Chiclayo  
 JOSE VIDAL VALLADOLID SALAZAR  
 FISICO MEDICO  
 LIC. IND. IPEN/OTAN N° 0927-27

  
 Clinicas **auna**/Sede Chiclayo  
 LUIS ENRIQUE QUISPE AYME  
 Dosimetrista en Teleterapia.  
 Lic. IPEN : 0268-20

  
 Clinicas **auna** | Sede Chiclayo  
 RICARDO PALMA ESPARZA  
 FISICO MEDICO  
 LIC. IND. IPEN/OTAN N° 8970-22

(iv) Recombination correction (two-voltage method)

Polarizing voltages:  $V_1$  (normal) = -300 V  $V_2$  (reduced) = -150 V  
Readings at each V:  $M_1$  = 12.788  $M_2$  = 12.746  
Voltage ratio  $V_1 / V_2$  = 2.0000 Read. ratio  $M_1 / M_2$  = 1.0033  
Beam type:   
 $a_0$  = 2.3370  $a_1$  = -3.6360  $a_2$  = 2.2990

1.0032<sup>d</sup>

Corrected dosimeter reading at the voltage  $V_1$ :

1.3075E-01

## 5. Absorbed dose to water at the reference depth, $z_{ref}$

Beam quality correction factor for user quality Q:

If  $Q_0$  is  $^{60}\text{Co}$  Table 18 gives

0.913

If  $Q_0$  is electron beam Table 19 gives

Use derived from

If is derived from series of electron beam calibrations

Calibration laboratory:

Date:

9.9134E-03 Gy / MU

## 6. Absorbed dose rate to water at the depth of dose maximum, $z_{max}$

Depth of dose maximum:  $z_{max}$  = 2.40 g cm<sup>-2</sup>

Percentage depth-dose at  $z_{ref}$  for a 10X10 cm x cm field size:

$PDD(z_{ref} = 2.8 \text{ g cm}^{-2})$  = 99.24 %

Absorbed-dose calibration of monitor at  $z_{max}$ :

9.9893E-03 Gy / MU

Notes:

300	-150	-300	Voltaje nC
-12.783	12.738	12.776	
-12.8	12.748	12.792	
-12.8	12.751	12.797	

-12.794 12.746 12.788

  
Clínicas **auna** | Sede Chiclayo  
JOSÉ VIDAL VALLADOÑO SALAZAR  
FÍSICO MÉDICO  
LIC. IND. IPEN/OTAN N° 0927-27

  
Clínicas **auna**/Sede Chiclayo  
LUIS ENRIQUE QUISPE AYME  
Dosimetrista en Teleterapia.  
Lic. IPEN : 0268-20

  
Clínicas **auna** | Sede Chiclayo  
RICARDO PALMA ESPARZA  
FÍSICO MÉDICO  
LIC. IND. IPEN/OTAN N° 8970-22

<sup>a</sup> Note that if  $Q_0$  is  $^{60}\text{Co}$ , is denoted  $N_{D,w}$

<sup>b</sup> All readings should be checked for leakage and corrected if necessary

<sup>c</sup>  $M$  in the denominator of  $k_{pol}$  denotes reading at the user polarity. Preferably, each reading in the equation should be the average of the ratios of  $M$  (or  $M_+$  or  $M_-$ ) to the reading of an external monitor,  $M_{em}$ .


<sup>d</sup> Check that

0.003

0.003

  
Clínicas **auna** | Sede Chiclayo  
JOSÉ VIDAL VALLADO SALAZAR  
FÍSICO MÉDICO  
LIC. IND. IPEN/OTAN N° 0922-22

  
Clínicas **auna**/Sede Chiclayo  
LUIS ENRIQUE QUISPE AYME  
Dosimetrísta en Teleterapia.  
Lic. IPEN : 0268-20

  
Clínicas **auna** | Sede Chiclayo  
RICARDO PALMA ESPARZA  
FÍSICO MÉDICO  
LIC. IND. IPEN/OTAN N° 8970-22